


# Dynamic Stress Rheometer (DSR200)



The DSR is a stress-controlled rheometer suitable for steady, dynamic and transient testing of both melts and fluids. Stresses are applied through a low-inertia AC eddy current motor attached to the upper geometry. The resulting angular velocities are measured by means of an optical encoder.

## Geometries

A broad variety of measurement geometries is available in the laboratory: plate-plate ( $\varnothing 25 / \varnothing 40$ , plastic and aluminium); cone-plate ( $\varnothing 25-0.1 / \varnothing 40-0.04$ , plastic and aluminium); couette cell (aluminium); disposable plates ( $\varnothing 25$ )

## Temperature control

Three different temperature control units are available: a recirculating fluid bath, a Peltier system (allowing temperature ramps up to  $50^{\circ}\text{C}/\text{min}$ ) and an oven that consist of electrically heated plates with an optional  $\text{N}_2$ -purge for temperatures up to  $350^{\circ}\text{C}$ .

## Instrument specifications

The instrument has a torque range of 0.001 up to 20 mNm with a torque resolution of 0.0001mNm at 0.1 mNm. The instrument can achieve a rotational rate up to 954 rpm and has an angular resolution of  $4.2 \mu\text{rad}$ . Its dynamic frequency range is between  $10^{-5}$  and 500 rad/s. The machine can be used as well to perform parallel superposition measurements.

Orchestrator software (Rheometrics) is used as an interface between user and instrument.



## Typical measurements

A broad spectrum of measurements can be performed with this instrument: e.g. the determination of yield stresses, monitoring the kinetics of crystallization, curing, ... using different time-temperature profiles, the determination of the morphology of blends by means of dynamic testing after different preshear histories, ...

